




MEDICAL EQUIPMENT
Haemodynamic monitoring



most-care^{Up}
up to the beat



Value Life

MostCare^{UP} is the only monitor able to follow, in real time and from beat to beat, even the slightest haemodynamic variations in the patient.

The patented algorithm based on the **PRAM method** (Pressure Recording Analytical Method) evaluates the cardiac output and many other haemodynamic parameters without any prior calibration.

An immediate, customisable interface can display a broad set of information regarding the preload, afterload, cardiac contractility and efficiency, which have become vitally important elements in optimising the treatment of high-risk patients and in defining the best haemodynamic settings for patients with alterations in their cardiovascular systems.



Advantages

Simple

- No calibration needed.
- Intuitive, customisable interface.
- No change in protocols in use.

Quick

- Constant monitoring with immediate results.
- Rapid connection and set up.

Versatile

- Any peripheral or femoral artery.
- Applicable to the widest range of patients.
- Easily transferred from one patient to another.

Innovative

- Exclusive variables (CCE, dirotic pressure, Ea).
- Patented dynamic filter to guarantee the quality of the pressure signal.
- Modern connectivity and data transfer systems.

Reliable

- Patented and validated algorithm.
- Immediate response to even the smallest haemodynamic variations.
- A wide range of clinical papers available.

Convenient

- Can be used on numerous patients without disposable nor added elements.
- *On Demand* system adaptable to all uses.

MostCare^{UP} is a reliable and efficient system that adapts to a wide range of patient types and clinical conditions.

Thanks to its rapid set-up, information can be obtained in real time and saved, reviewed and transferred for subsequent analysis.




The *Endless* and *On Demand* versions allow the operator to choose the method of use which best suits the specific needs, thereby guaranteeing an effective control on costs.

PRAM method

A patented algorithm

PRAM (Pressure Recording Analytical Method) is an innovative method to analyse the pressure wave used in MostCare^{Up}¹. It allows for constant and sensitive monitoring in real time of the slightest haemodynamic variations because it is based, heart beat by heart beat, only on the morphology of the arterial pressure wave.

- Sampling at 1000 Hz
- Beat by beat analysis of the wave form
- Does not depend on pre-estimates
- No external calibration required




Atrial fibrillation Aortic counterpulsation Fluid responsiveness Aortic Insufficiency

Vasodilatation Vasoconstriction

Each patient is unique and his haemodynamic condition can evolve rapidly. The shape of the arterial pressure wave is the result of a complex balance which depends on both the coupling of the cardiac function with the vascular system and their interaction with the respiratory system.


The precise analysis of the shape of the wave obviates the need for calibration and preestimated data about the patient. It also identifies the dicrotic pressure and the Z(t) impedance of the cardiovascular system, even in cases of unusual pressure wave forms.



¹ Romano SM, Pistolesi M, Crit Care Med, 2002


Haemodynamic variables

Pressure




Systolic, diastolic, mean and pulse pressure (PP) are measured with every heartbeat.

Dicrotic pressure




The value of the dicrotic pressure, gauged with precision at 1000Hz, provides information about the vascular condition and the ventricle-arterial coupling.

Cardiac output




The stroke volume (SV) is measured beat-by-beat and allows for the cardiac output (CO) to be calculated.


Derived variables



Systemic vascular resistance (SVR), cardiac power output (CPO) and oxygen delivery (DO₂) are examples of the derived variables provided by MostCare^{Up}.



CCE




Cardiac cycle efficiency (CCE) is an exclusive variable which describes haemodynamic performance in terms of energy expenditure in the patient being monitored².

dP/dt_{max}




The maximum pressure variation compared to time (dP/dt_{max}) is linked to the heart's contractility and also to the condition of the vascular system.

Dynamic variables



Pulse pressure variation (PPV) and stroke volume variation (SVV) during the respiratory cycle can be viewed simultaneously.

Dynamic filter




The shape of the pressure curve can be affected by resonance phenomena. The exclusive, dynamic filter in MostCare^{Up} has been designed to automatically optimise the quality of the wave and to reduce these phenomena³.

² Romano SM, Int J Cardiol, 2012
³ Romagnoli S et al., Crit Care, 2014

Markers and trends Do & check function

MostCare^{Up} can display trends for many haemodynamic variables simultaneously. It is also possible to insert personalised markers during specific events (e.g. start treatment). The **do&check function** was designed specifically to help the clinic when monitoring haemodynamic variations following specific treatments (e.g. fluid challenge).



Connectivity and data management

MostCare^{Up} supports the most advanced communication and data transmission standards. The patient's parameters and freezer-frames can be saved to the machine's memory or exported via the USB port. The data can also be transferred to the hospital's platform using the HL7 protocol. The image on the display can be shared for monitoring or educational purposes via HDMI.



Flexibility of use *Endless-On Demand*

MostCare^{Up} guarantees the maximum flexibility and cost efficiency thanks to the various ways in which it can be used. The *On Demand* version of the monitor can be activated for a single use or for periods of time to meet specific application needs. The *Endless* version allows an unlimited use of the system without additional cost.

Applications

The PRAM method requires no external calibration or anthropometric normalisation.

MostCare^{Up} can therefore be easily used on any patient who requires constant or occasional haemodynamic monitoring. More specifically, during haemodynamic instability or in the presence of acute clinical variations in high risk patients.

Goal directed therapy

Perioperative

Fluid optimisation in high risk surgery patients has significantly reduced postoperative complications, length of hospital stay and improved the outcome of the patients, thus resulting in substantially lower costs.

⁴ Pearse R *et al.*, Crit Care, 2005

⁵ Lopes MR *et al.*, Crit Care, 2007

⁶ Vincent JL *et al.*, Crit Care, 2015

Intensive care and critical patients

Thanks to beat by beat analysis, the PRAM method is able to reliably recognise and monitor the haemodynamic changes resulting from the administration of vasoactive drugs and fluids, in real time and even in septic or trauma patients.

⁷ Vincent JL *et al.*, Crit Care, 2011

⁸ Franchi F *et al.*, BJA, 2011

⁹ Guarracino F *et al.*, Crit Care, 2014

¹⁰ Donati A *et al.*, J Crit Care, 2014

Evaluation of the ventricular function

Critical patient

The echocardiography is a gold standard technique to evaluate ventricular function. Some of the variables provided by MostCare^{Up} (dP/dt_{max} and CCE) supply a constant stream of information about the cardiac function in the critical patient.

¹¹ Scolletta S *et al.*, Intensive Care Med, 2013

Cardiac insufficiency

The beat by beat monitoring of haemodynamic variables like diastolic pressure, dP/dt_{max} and CCE guarantees a rapid and immediate evaluation of any clinical variations in the patient so that immediate action can be taken.

¹² Giglioli C *et al.*, Eur J Heart Fail, 2011

¹³ Pavoni V *et al.*, J Anesth Clin Res, 2012

¹⁴ Barile L *et al.*, J. Cardiothorac Vasc Anesth, 2013



Specific applications

Paediatric patient

¹⁵ Calamandrei M et al., *Pediatr Crit Care Med*, 2008

¹⁶ Ricci Z et al., *Crit Care*, 2014

¹⁷ Garisto C et al., *Paediatr Anaesth*, 2014

Ventilation

¹⁸ McBride WT et al., *J Cardiothorac Vasc Anesth*, 2012

Obese patient

¹⁹ Balderi T et al., *Obes Surg*, 2008

Aortic counterpulsation

²⁰ Zangrillo A et al., *J Cardiothorac Vasc Anesth*, 2010

²¹ Gelsomino Set al., *Eur J Cardiothorac Surg*, 2012

²² Onorati F et al., *J Thorac Cardiovasc Surg*, 2012


Interventional cardiology

²³ Romagnoli S et al., *J Cardiothorac Vasc Anesth*, 2010

²⁴ Giglioli C et al., *World J Cardiovasc Dis*, 2013

Hypothermia therapy

²⁵ Lazzeri C et al., *Acute Card Care*, 2014



Variables of MostCare^{Up}

Haemodynamic variables	Formulas	Physiological range ***	Units
Pressures			
Sys	Systolic pressure		mmHg
Dia	Diastolic pressure		mmHg
MAP	Mean arterial pressure		mmHg
Dic	Dicotic pressure	70 - 105	mmHg
PP	Pulse pressure	$P_{sys} - P_{dia}$	30 - 50 mmHg
MAP-Dic	Mean and dicotic pressure difference	MAP-Dic	-10 - +10 mmHg
CVP*	Central venous pressure		mmHg
Cardiac output			
SV	Stroke volume		60 - 100 mL
SVI	Stroke volume index		35 - 45 mL/m ²
SV _{kg}	Weighted stroke volume	SV/weight	mL/kg
CO	Cardiac output		4.0 - 8.0 L/min
CI	Cardiac output index		2.6 - 3.8 L (min · m ²)
SVR	Systemic vascular resistance	$(MAP - CVP) / CO \cdot 80$	800 - 1400 dyne · sec/cm ⁵
SVRI	Systemic vascular resistance index	$(MAP - CVP) / CI \cdot 80$	1600 - 2400 dyne · sec · m ² /cm ⁵
Oxygen delivery			
SpO ₂ *	Arterial oxygen saturation		96 - 100 %
DO ₂ *	Oxygen delivery	$DO_2 = CO \cdot CaO_2$ con $CaO_2 = Hb \cdot 1,34 \cdot SaO_2$	900 - 1000 mL/min
DO ₂ I*	Oxygen delivery index	$DO_2I = DO_2 / BSA$	500 - 600 mL/min/m ²
Efficiency and cardiac function			
dP/dt _{max}	Maximal slope of the systolic upstroke		0.9 - 1.3 mmHg/msec
CCE	Cardiac cycle efficiency		-0.2 - 0.3 units
CPO	Cardiac power	$MAP \cdot CO / 451$	0.80 - 1.20 W
CPI	Cardiac power index	$MAP \cdot CI / 451$	0.50 - 0.70 W/m ²
Vascular function			
Ea	Arterial elastance	Dic/SV	1.10 - 1.40 mmHg/mL
PPV/SVV	Dynamic elastance	PPV/SVV	units
Z _{tot}	Cardiovascular impedance		mmHg · sec/mL
Dynamic variables			
PPV	Pulse pressure variation		< 15** %
SVV	Stroke volume variation		< 15** %
SPV	Systolic pressure variation		%
DPV	Dicotic pressure variation		%
Other specific variables			
PR	Pulse rate		1/min
Dia _{pk}	Diastolic peak		mmHg

*When added probes are used. DO₂ and DO₂I calculated with fixed Hb value.

** Approximate values reported in the literature in the patient receiving controlled mechanical ventilation.

*** Normal values in the adult patient. The values depend on the patient in relation to the clinical conditions.

BSA = body surface area, calculated by the standard formulas of DuBois & DuBois, using the values of weight and height.

Bibliography

- Romano SM, Pistolesi M. Assessment of cardiac output from systemic arterial pressure in humans. *Crit Care Med* 2002 Aug; 30(8): 1834-41.
- Romano SM. Cardiac cycle efficiency: A new parameter able to fully evaluate the dynamic interplay of the cardiovascular system. *Int J Cardiol* 2012 Mar; 155(2): 326-7.
- Romagnoli S, Ricci Z, Quattrone D, Tofani L, Tujjar O, Villa G, Romano SM, De Gaudio A. Accuracy of invasive arterial pressure monitoring in cardiovascular patients: an observational study. *Crit Care* 2014 Nov 30; 18(6): 644.
- Pearse R, Dawson D, Fawcett J, Rhodes A, Grounds RM, Bennett ED. Early goal-directed therapy after major surgery reduces complications and duration of hospital stay. A randomised, controlled trial. *Crit Care* 2005; 9(6): R687-93.
- Lopes MR, Oliveira MA, Pereira VO, Lemos IP, Auler JO Jr, Michard F. Goal-directed fluid management based on pulse pressure variation monitoring during high-risk surgery: a pilot randomized controlled trial. *Crit Care* 2007; 11(5): R100.
- Vincent JL, Pelosi P, Pearse R, Payen D, Perel A, Hoefl A, Romagnoli S, Ranieri VM, Ichai C, Forget P, Della Rocca G, Rhodes A. Perioperative cardiovascular monitoring of high-risk patients: An international consensus. *Crit Care* 2015 May 8; 19: 224.
- Vincent JL, Rhodes A, Perel A, Martin GS, Della Rocca G, Vallet B, Pinsky MR, Hofer CK, Teboul JL, de Boode WP, Scolletta S, Vieillard-Baron A, De Backer D, Walley KR, Maggiorini M, Singer M. Clinical review: Update on hemodynamic monitoring. A consensus of 16. *Crit Care* 2011 Aug; 15(4): 229.
- Franchi F, Silvestri R, Cubattoli L, Taccone FS, Donadello K, Romano SM, Giomarelli P, McBride WT, Scolletta S. Comparison between an uncalibrated pulse contour method and thermodilution technique for cardiac output estimation in septic patients. *Br J Anaesth* 2011; 107(2): 202-8.
- Guarracino F, Ferro B, Forfori F, Bertini P, Magliacane L, Pinsky MR. Jugular vein distensibility predicts fluid responsiveness in septic patients. *Crit Care* 2014 Dec 5; 18(6): 647.
- Donati A, Carsetti A, Tondi S, Scorcella C, Domizi R, Damiani E, Gabbanelli V, Münch C, Adrario E, Pelaia P, Cecconi M. Thermodilution vs pressure recording analytical method in hemodynamic stabilized patients. *J Crit Care* 2014 Feb; 29(2): 260-4.
- Scolletta S, Bodson L, Donadello K, Taccone FS, Devigili A, Vincent JL, De Backer D. Assessment of left ventricular function by pulse wave analysis in critically ill patients. *Intensive Care Med* 2013; 39: 1025-33.
- Giglioli C, Landi D, Cecchi E, Chiostrini M, Gensini GF, Valente S, Ciaccheri M, Castelli G, Romano SM. Effects of ULTRAFiltration vs DiureticS on clinical neurohormonal and hemodynamic variables in patients with deCompensated heart failure: the ULTRADISCO study. *Eur J Heart Fail* 2011; 13(3): 337-46.
- Pavoni V, Romagnoli S, Batignani G, Gianesello L, Horton A, Romano SM. Unsuspected Heart Failure: Usefulness of a Minimally Invasive Hemodynamic Monitoring System. *J Anesth Clin Res* 2012; 3(8).
- Barile L, Landoni G, Pieri M, Ruggeri L, Maj G, Nigro Neto C, Pasin L, Cabrini L, Zangrillo A. Cardiac Index Assessment by the Pressure Recording Analytic Method in Critically Ill Unstable Patients After Cardiac Surgery. *J Cardiothorac Vasc Anesth* 2013 Dec; 27(6): 1108-13.
- Calamandrei M, Mirabile L, Musichetta S, Gensini GF, De Simone L, Romano SM. Assessment of cardiac output in children: a comparison between the pressure recording analytical method (PRAM) and the Doppler echocardiography method - a pilot study. *Pediatr Crit Care Med* 2008; 9: 310-2.
- Ricci Z, Haiberger R, Pezzella C, Garisto C, Favia I, Cogo P. Furosemide versus ethacrynic acid in pediatric patients undergoing cardiac surgery: a randomized controlled trial. *Crit Care*. 2015 Jan 7; 19(1): 2.
- Garisto C, Favia I, Ricci Z, Romagnoli S, Haiberger R, Polito A, Cogo P. Pressure recording analytical method and bioreactance for stroke volume index monitoring during pediatric cardiac surgery. *Paediatr Anaesth* 2015 Feb; 25(2): 143-9.
- McBride WT, Ranaldi G, Dougherty MJ, Siciliano T, Trethowan B, Elliott P, Rice C, Scolletta S, Giomarelli P, Romano SM, Linton DM. The hemodynamic and respiratory effects of continuous negative and control-mode cuirass ventilation in recently extubated cardiac surgery patients: part 2. *J Cardiothorac Vasc Anesth* 2012; 26: 873-7.
- Balderi T, Forfori F, Marra V, Di Salvo C, Dorigo M, Anselmino M, Romano SM, Giunta F. Continuous Hemodynamic Monitoring During Laparoscopic Gastric Bypass in Superobese Patients by Pressure Recording Analytical Method. *Obes Surg* 2008 Aug; 18(8): 1007-14.
- Zangrillo A, Maj G, Monaco F, Scandroglio AM, Nuzzi M, Plumari V, Virzo I, Bignami E, Casiraghi G, Landoni G. Cardiac Index Validation Using the Pressure Recording Analytic Method in Unstable Patients. *J Cardiothorac Vasc Anesth* 2010; 24: 265-9.
- Gelsomino S, Renzulli A, Rubino AS, Romano SM, Lucà F, Valente S, Gensini GF, Lorusso R. Effects of 1:1, 1:2 or 1:3 intra-aortic balloon counterpulsation/heart support on coronary haemodynamics and cardiac contractile efficiency in an animal model of myocardial ischaemia/reperfusion. *Eur J Cardiothorac Surg* 2012; 42: 325-32.
- Onorati F, Santini F, Amocelli E, Campanella F, Chiominto B, Faggian G, Mazzucco A. How should I wean my next intra-aortic balloon pump? Differences between progressive volume weaning and rate weaning. *J Thorac Cardiovasc Surg* 2012; 145: 1214-21.
- Romagnoli S, Romano SM, Bevilacqua S, Lazzeri C, Santoro G, Ciappi F, Gelsomino S, Dini D. Pulse contour cardiac output monitoring during a complicated percutaneous aortic valve replacement. *J Cardiothorac Vasc Anesth* 2010 Apr; 24(2): 303-5.
- Giglioli C, Tujjar O, Cecchi E, Landi D, Chiostrini M, Valente S, Baldereschi GJ, Meucci F, Gensini GF, Romano SM. Hemodynamic changes acutely determined by primary PCI in STEMI patients evaluated with a minimally invasive method. *World J Cardiovasc Dis* 2013; 3: 69-72.
- Lazzeri C, Sori A, Bernardo P, Chiostrini M, Tommasi E, Zucchini M, Romano SM, Gensini GF, Valente S. Cardiovascular effects of mild hypothermia in post-cardiac arrest patients by beat-to-beat monitoring. A single centre pilot study. *Acute Card Care*. 2014 Jun; 16(2): 67-73

Product codes

Monitors	
VMB08MC0202E0V	MostCare ^{Up} monitor – On-Demand version - Standard
VMB08MC0255E0V	MostCare ^{Up} monitor – Endless version - Standard
VMB08MC0202E1V	MostCare ^{Up} monitor – On-Demand version - Turkey
VMB08MC0255E1V	MostCare ^{Up} monitor – Endless version - Turkey
VMB08MC0202E2V	MostCare ^{Up} monitor – On-Demand version - UK
VMB08MC0255E2V	MostCare ^{Up} monitor – Endless version - UK
VMB08MC0202E3V	MostCare ^{Up} monitor – On-Demand version - AUS/NZ/VN
VMB08MC0255E3V	MostCare ^{Up} monitor – Endless version - AUS/NZ/VN

Cards	
VMB04MCU1D02	MostCare ^{Up} card – 24 hours (1 day)
VMB04MCU1M02	MostCare ^{Up} card – 30 days (1 month)
VMB04MCU1Y02	MostCare ^{Up} card – 365 days (1 year)
VMB04MCU01U02	MostCare ^{Up} card – 1 use
VMB04MCU10U02	MostCare ^{Up} card – 10 uses
VMB04MCU50U02	MostCare ^{Up} card – 50 uses

Cables	
Direct Cables	
VMC03MUBBR	MostCare ^{Up} direct BP cable – B.BRAUN type
VMC03MUBD	MostCare ^{Up} direct BP cable – BD type
VMC03MUBIO	MostCare ^{Up} direct BP cable – BIOSENSOR/UTAH type
VMC03MUDPT	MostCare ^{Up} direct BP cable – CODAN BDPT type
VMC03MUEDW	MostCare ^{Up} direct BP cable – EDWARDS type
VMC03MUMED	MostCare ^{Up} direct BP cable – MEDEX LOGICAL type
VMC03MUMED2	MostCare ^{Up} direct BP cable – MEDEX TRANSTAR type
Y cables	
VMC03MUYBBR	MostCare ^{Up} Y BP cable – B.BRAUN type
VMC03MUYBD	MostCare ^{Up} Y BP cable – BD type
VMC03MUYBIO	MostCare ^{Up} Y BP cable – BIOSENSOR/UTAH type
VMC03MUYDPT	MostCare ^{Up} Y BP cable – CODAN BDPT type
VMC03MUYEDW	MostCare ^{Up} Y BP cable – EDWARDS type
VMC03MUYMED2	MostCare ^{Up} Y BP cable – MEDEX TRANSTAR type
Analog Cables	
VTCMUHEMO018	Input/Output Analog Cable Jack 3,5
VTCMUHEMO024	Input/Output Analog Cable DS
VTCMUHEMO031	Input/Output Analog Cable 7 Pin
VTCMUHEMO033	Input/Output Analog Cable PDM
VTCMUHEMO035	Input/Output Analog Cable Jack 6

 **CRITICAL CARE**

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only
and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE
Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81
Export customer service: +33 (0)1.39.92.64.15
Fax.: +33 (0)1.39.92.64.44 • www.vygon.com

 @VygonGroup  Vygon  VygonGroup  VygonGroup



www.vygon.com