




EMERGENCY
Difficult Intubation




bO₂ugie **Boussignac**
ET tube introducer



Value Life

● Easier endotracheal intubation

bO₂ugie **Boussignac** is a **single-use introducer** which makes endotracheal intubation easier in patients with difficult intubation.




Endotracheal tube introducer allowing oxygen delivery

If oxygen delivery is needed, attach connector at the proximal end of the bougie.



Side holes for decompression

In case of **laryngospasm**, during oxygen delivery, oxygen can escape through the side holes preventing lung damage (e.g. volutrauma/barotrauma).



Depth markings for accurate insertion



Code	Size	Colour	Length	Endotracheal tube inner diameter	Sales Unit
5576.01	15Fr (Adult)	Green	70cm	6mm to 10mm	10 units/box
5576.02	10Fr (Children)	Blue	60cm	4mm to 5.5mm	
5576.03	6Fr (Newborn)	Pink	50cm	2.5mm to 3.5mm	

All products supplied with connectors.



ANAESTHESIA & EMERGENCY

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE
 Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81
 Export customer service: +33 (0)1.39.92.64.15
 Fax.: +33 (0)1.39.92.64.44 • www.vygon.com



www.vygon.com