



VASCULAR ACCESS

Silicone catheter for long-term use



lifecath

Long-term venous catheter
with cuff



Value Life

● Silicone catheter with subcutaneous cuff for long-term venous access.

Radiopaque silicone catheter for long-term venous access

(parenteral nutrition, administration of medication, antibiotics, chemotherapy, etc.).

- Inserted by venous cut-down (catheter only) or percutaneous puncture (kit with peel-apart Desilet)
- Subcutaneous cuff to fix the catheter securely to subcutaneous tissues
- Reinforced extra-vascular section with clamping site and clamp
- Fixed proximal hub therefore tunnelling is carried out before venipuncture
- Tunnelling needle supplied with all kits

Silicone:
Haemocompatibility
Biostability
Superior flexibility

lifecath



Ordering

• lifecath

Code		Catheter					Peel-apart desilet
Catheter only	Catheter With introducer	Ø Fr	Flow rate ml/min	Priming volume ml	Total length cm	Useful length cm	Fr
2191.27		2.7	1.0	0.25	75	50	
2191.42	2191.425	4.2	5.0	0.35	75	50	5
2191.50	2191.506	5	15.0	0.70	75	50	6
2191.66	2191.667	6.6	27.0	1.10	90	55	7
2193.96	2193.960	9.6	60.0	1.65	90	55	10

repair kit (!)	
Code	
2194.27	
2194.42	
2194.50	
2194.66	
2195.96	

• Double-lumen lifecath

Code		Catheter				
Catheter only	Catheter With introducer	Ø Fr	Flow rate ml/min	Priming volume ml	Total length cm	Useful length cm
2293.070	2293.70	7	2 - 15	0.70 – 1.10	90	55
2293.095	2293.95	9	37 - 37	1.44 – 1.44	90	55
2293.110	2293.11	11	25 - 80	1.30 – 2.20	90	55
2293.140	2293.14	14	80 - 80	2.20 – 2.20	90	55

repair kit (!)	
Code	
For extension tube*	Complete Kit**
2195.06/10	2295.07
2195.13	2295.95
2195.10/16	2295.11
2195.16	2295.14

*For repairing one extension tube, proximal to the fixation hub

**For repairing the external part of the catheter (distal to the fixation hub), this model includes 2 extension tubes

- Flowrate and priming volume are calculated for the complete catheter before the length is adjusted by the physician.

(!) When ordering the repair **kit**, the silicone glue must be ordered separately under code number 2197.000.

Minimum ordering quantity: **1 piece**

 INTRAVASCULAR THERAPIES

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE
Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81
Export customer service: +33 (0)1.39.92.64.15
Fax: +33 (0)1.39.92.64.44 • www.vygon.com


www.vygon.com