





The **high security** connection







Value Life

Preparation devices for medications/solutions

Bag


Brand name	Code	Product	Description	Indication	Quantity box/case
qimospike B	7218.01		qimospike B is an IV fluid bag access spike fitted with qimo ♀.	Aspiration or addition of medication/ fluid from/into an IV fluid bag.	50/600
qimospike B valve	7218.11		qimospike B valve is an IV fluid bag access spike fitted with qimo ♀ and a non-return valve.	Aspiration only from an IV fluid bag, thus avoiding the accidental contamination of the IV solution.	25/450

Vial

Brand name	Code	Product	Description	Indication	Quantity box/case
qimocap ø14	7219.01 7219.010*		qimocap is a vial access cap fitted with qimo ♀. It has been specifically designed to create a permanent connection with the medication vial.	Withdrawal and reconstitution of small volumes of medication in a completely closed environment.	50/900 4/72
qimocap ø20	7219.02 7219.020*				25/450 4/48
qimospike V	7217.04 7217.040*		qimospike V is a vented spike for vials fitted with qimo ♀.	Repeated withdrawals from every type of vial.	25/300 4/48
qimoharpoon	7213.02		qimoharpoon is a closed system transfer device that has an air-intake tube connected to a 50ml balloon filled with sterile air which ensures pressure equalization.	Withdrawal and reconstitution of high-risk medications.	10/80
qimoharpoon S	7213.01				

*Multipack by 25

Syringe/Extension tube

Brand name	Code	Product	Description	Indication	Quantity box/case
qimo ♂	7210.91 7210.910*		qimo ♂ is a disinfectable closed needleless safety connector. It is a safety captive adaptor for male Luer-lock hubs. Once attached, it cannot be disconnected. It connects only to qimo ♀.	To secure all the male Luer-lock fittings.	50/900 4/72
qimotransfer	7220.01		Female Luer-lock to female Luer-lock connector, with qimo ♀ on both ends of the connector.	To transfer medications from syringe to syringe.	50/600



*Multipack by 25

Elastomeric pump - 2 packs available

	Brand name	Code	Product	Description	Indication	Quantity box/case
Option 1	qimoprime, qimo ♂ & qimo ♀	7210.04		Set with qimoprime, qimo ♂ and qimo ♀ in the same packaging.	Can be used to secure the filling port of the elastomeric pump and the distal end of its tubing, but other applications are possible.	50/600
Option 2	qimo ♀	7210.02		qimo ♀ is a multi-purpose, closed, needleless and disinfectable safety connector, supplied in individual packaging.	qimo ♀ can be used to connect the tubing of the distal end of the elastomeric pump secured with qimo ♂, to access an implantable port or any other patient vascular device.	50/1200
	qimoprime & qimo ♂	7210.01		Set supplied with qimo ♂ and qimoprime in the same packaging.	Can be used to secure distal end of elastomeric pump tubing. Need to be connected to a qimo ♀	50/600

Transfer devices for medications/solutions







Bag

Brand name	Code	Product	Description	Indication	Quantity box/case
qimconnect	7221.01 7221.020*		IV fluid bag extension set with a lateral injection site fitted with qimo ♀. qimo ♂ secures the distal end of the extension tube.	Dilution and administration of all kind of medications.	25/150 10/60
qimconnect filter	7221.011 7221.0110*		IV fluid bag extension set with a 0.2µm filter and a lateral injection site fitted with qimo ♀. qimo ♂ secures the distal end of the extension tube.	Dilution and administration of solutions which need to be filtered (ex: Paclitaxel).	25/150 10/60




*Multipack by 10 (without qimo ♂)

Infusion devices




Generic products

Brand name	Code	Product	Description	Indication	Quantity box/case
qimopump	7227.101 7227.201 7227.401		qimopump is a safe adaptor for pump infusion sets. This device is fitted with 1, 2 or 4 qimo ♀ connectors, a spike with air-venting port (with a cap) and a universal tube to connect to pump dedicated sets.	To secure the dedicated pump infusion sets and connect the bag containing medication in a safe manner.	25/150
qimo IVset with integral qimooctopus	7229.101 7229.201 7229.401		Intravenous infusion device with spike and drip chamber. It is fitted with 1, 2 or 4 qimo ♀ and qimo ♂ to attach to the male Luer-lock port at the distal end of qimo IVset. qimoprime allows the priming of the line once qimo ♂ has been attached.	To infuse safely medications by gravity or with qimono infusion pump (QIP).	25/150
qimo IVset without qimooctopus	7229.102 7229.202 7229.402				
qimconnect ♀	7221.02		IV fluid bag extension tube with a lateral injection site fitted with qimo ♀ and an air-venting port (with a cap).	Preparation and administration of non toxic drugs (premedication such as antiemetics, corticoids, etc.)	25/150
qimoprime & qimo ♂	7210.01		Set supplied with qimo ♂ and qimoprime in the same packaging.	To secure both the male Luer-lock and the female Luer-lock end of infusion devices.	50/600
qimooctopus	7231.221		qimooctopus is a safe double-lumen IV extension tube allowing continuous infusion, intermittent injection and blood sampling. Its lateral line is fitted with qimo ♀. qimo ♂ has to be attached to the distal end of qimooctopus. qimoprime allows the priming of the line once qimo ♂ has been attached.	Safe emergency access. Easy direct injection with a syringe close to the patient.	50/300

Infusion sets and Qimono Infusion Pump (QIP)

Brand name	Code	Product	Description	Indication	Quantity box/case
Infusion set for gravity or peristaltic pump					
qimo IVset with integral qimo octopus	7229.101 7229.201 7229.401		Intravenous infusion device with spike and drip chamber. It is fitted with 1, 2 or 4 qimo ♀ and qimo ♂ to attach to the male Luer-lock port at the distal end of qimo IVset .	To infuse safely medications by gravity or with qimono infusion pump (QIP).	25/150
qimo IVset without qimo octopus	7229.102 7229.202 7229.402		qimo prime allows the priming of the line once qimo ♂ has been attached.		
Accessory					
qimono infusion pump (QIP)	0VEPQIM FR 0VEPQIMG02 EN PLUG 0VEPQIMC03 NL 0VEPQIMC04 DE 0VEPQIMC06 ES 0VEPQIMC07 PT 0VEPQIMC08 IT 0VEPQIMC10 PL		Infusion Pump : qimono infusion pump (QIP) is safe, simple to use, precise, silent. 2 set up modes available : Volume & Time of infusion or Flow rate & Volume to infuse. It displays a large coloured screen, visual + audible critical alarms for nurse monitoring & patient safety. It is easy to clean and offers a quick access to the Lithium battery. 2 years warranty (excluding battery).	Qimono Infusion Pump is dedicated to administration of cytotoxic drugs and perfusion in day care hospital. To be used ONLY in conjunction with qimo IVset listed above	1 unit/box

Patient's vascular access

Brand name	Code	Product	Description	Indication	Quantity box/case																																															
qimo perfusafe	7242.720 7242.725 7242.920 7242.925		qimo perfusafe is a safety infusion set with Huber needle fitted with a qimo ♀ to access implantable ports.	Safety device (protected needle and closed system connector) to puncture an implantable port.	10/90																																															
qimo ♀	7210.02		qimo ♀ is a multi-purpose, closed, needleless and disinfectable safety connector.	To connect to any other vascular access, Huber needles, CVC...safety and with a closed system	50/1200																																															
qimo flow+	<table border="1"> <thead> <tr> <th>Gauge</th> <th>Ø</th> <th>Longueur</th> <th>Code</th> </tr> </thead> <tbody> <tr> <td rowspan="5">22G</td> <td rowspan="5">0,7</td> <td>15mm</td> <td>741507*</td> </tr> <tr> <td>17mm</td> <td>741707</td> </tr> <tr> <td>20mm</td> <td>742007</td> </tr> <tr> <td>25mm</td> <td>742507</td> </tr> <tr> <td>30mm</td> <td>743007</td> </tr> <tr> <td rowspan="5">20G</td> <td rowspan="5">0,9</td> <td>15mm</td> <td>741509*</td> </tr> <tr> <td>17mm</td> <td>741709</td> </tr> <tr> <td>20mm</td> <td>742009</td> </tr> <tr> <td>25mm</td> <td>742509</td> </tr> <tr> <td>30mm</td> <td>743009</td> </tr> <tr> <td rowspan="5">19G</td> <td rowspan="5">1,1</td> <td>15mm</td> <td>741511*</td> </tr> <tr> <td>17mm</td> <td>741711</td> </tr> <tr> <td>20mm</td> <td>742011</td> </tr> <tr> <td>25mm</td> <td>742511</td> </tr> <tr> <td>30mm</td> <td>743011</td> </tr> <tr> <td></td> <td></td> <td>35mm</td> <td>743511</td> </tr> </tbody> </table>				Gauge	Ø	Longueur	Code	22G	0,7	15mm	741507*	17mm	741707	20mm	742007	25mm	742507	30mm	743007	20G	0,9	15mm	741509*	17mm	741709	20mm	742009	25mm	742509	30mm	743009	19G	1,1	15mm	741511*	17mm	741711	20mm	742011	25mm	742511	30mm	743011			35mm	743511		Safety huber needle with extension tube and safety adaptor: qimo ♀ (closed system) qimo flow+ is a single-handed activation needle allowing an automatic positive flush during removal	qimo flow+ safety Huber needles are indicated for the administration or withdrawal of fluids through implantable ports.	12/30
	Gauge	Ø	Longueur	Code																																																
	22G	0,7	15mm	741507*																																																
			17mm	741707																																																
			20mm	742007																																																
			25mm	742507																																																
			30mm	743007																																																
	20G	0,9	15mm	741509*																																																
			17mm	741709																																																
			20mm	742009																																																
			25mm	742509																																																
			30mm	743009																																																
	19G	1,1	15mm	741511*																																																
			17mm	741711																																																
20mm			742011																																																	
25mm			742511																																																	
30mm			743011																																																	
		35mm	743511																																																	

* Please check availability with your Vygon local representative

 INTRAVASCULAR THERAPIES

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only
and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE
Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81
Export customer service: +33 (0)1.39.92.64.15
Fax.: +33 (0)1.39.92.64.44 • www.vygon.com



www.vygon.com