



NUTRITION  
Enteral feeding

nutrifit®



Safe enteral feeding  
system for adults and  
paediatrics



Value Life

## What is ISO 80369 ?

It is an International Standard that aims to standardize all kinds of small-bore connectors for liquids and gases in healthcare applications.

The objective is to avoid misconnections between devices that have completely different functions (e.g. feeding tube connected to an IV line).

ISO 80369 consists of :

- Part 1: General requirements
- Part 2: Connectors for respiratory application
- **Part 3: Connectors for enteral application**
- Part 4: Connectors for urinary application
- Part 5: Connectors for limb cuff inflation application
- Part 6: Connectors for neuraxial application
- Part 7: Connectors for intravascular application (Luer)
- Part 20: Test methods



## What is ISO 80369-3 ?

**ISO 80369-3** is part 3 of ISO 80369, which is dedicated to enteral application.

This part 3 describes all the dimensions and drawings of the chosen ISO enteral connector.


Today, only one ISO enteral connector has been chosen for all patients, from Neonatology to Adult.

## What is « ENFit™ » ?

To simplify communication, enteral connectors compliant with ISO 80369-3 will be called « **ENFit™ connectors** ».

## What is nutri**fit**® ?

nutri**fit**® is the brand name of Vygon ENFit™ range.




## Is nutri**fit**® compliant with ISO 80369-3 ?

Yes, nutri**fit**® is compliant with ISO 80369-3.

## ● What is Vygon position concerning the ENFit™ system ?

Vygon recommends ENFit™ for paediatrics & adults.

A risk of under/over - delivery for small volumes exists especially in neonates.



## ● What does Vygon recommend for neonates ?

Vygon recommends nutri**safe**2® for neonatology. This safe enteral feeding system has been existing for 10 years. Safe, small and accurate, nutri**safe**2® is compliant with ISO 80369-1.



## ● nutri**fit**®: a safe enteral feeding system for adults and paediatrics

nutri**fit**® offers a full range of products:

- Pumps
- Sets for pumps
- Feeding tubes
- Syringes
- Extension sets
- Sampling straws
- Accessories



• **PVC tubes**

Length	Size					
	8 Fr	10 Fr	12 Fr	14 Fr	16 Fr	18 Fr
90cm	366.087	366.107	366.127	366.147	366.167	366.187
125cm	362.087	362.107	362.127	362.147	362.167	362.187

Without stylet

• **PUR tubes without stylet**

Length	Size					
	8 Fr	10 Fr	12 Fr	14 Fr	16 Fr	18 Fr
40cm	1364.087					
50cm	1361.087	1361.107				
90cm	1366.087	1366.107	1366.127			
125cm	1362.087	1362.107	1362.127	1362.147	1362.167	1362.187

• **PUR tubes with stylet**

Length	Size					
	8 Fr	10 Fr	12 Fr	14 Fr	16 Fr	18 Fr
90cm closed tip	1366.0871	1366.1071	1366.1271	1366.1471	1366.1671	1366.1871
90cm open tip	1366.0872	1366.1072	1366.1272			
125cm closed tip	1362.0871	1362.1071	1362.1271	1362.1471	1362.1671	1362.1871
125cm open tip	1362.0872	1362.1072	1362.1272			
125cm with weighted tip	1362.0874	1362.1074	1362.1274	1362.1474	1362.1674	1362.1874

• **Silicone tubes**

Length	Size					
	8 Fr	10 Fr	12 Fr	14 Fr	16 Fr	18 Fr
50cm	2331.087					
125cm	2332.087	2332.107	2332.127	2332.147	2332.167	2332.187

Without stylet

• **PVC dual flow gastric aspiration tube**

Length	Size					
	8 Fr	10 Fr	12 Fr	14 Fr	16 Fr	18 Fr
60cm	340.087					
90cm		340.107	340.127	340.147		
120cm					340.167	340.187

Without stylet

## nutrifit®: syringes

### • Transparent syringes

Code	Volume
1015.017	1ml
1015.027	2.5ml
1015.057	5ml
1015.107	10ml
1015.207	20ml
1015.607	60ml



## nutrifit®: accessories

### • Extension tubes

Length	Standard diameter (1.5 mm)
100cm	368.107
150cm	368.157

### • Extensions lines with clamp

Length	Standard diameter (1.5 mm)
100cm	368.1071
150cm	368.1571

### • Sampling straws

Length	Very small bore Ø 1.2x2.0mm
5cm	817.057
15cm	817.157

### • Other accessories

Product	Code
Y connector	368.47
Female cap only	829.07
3 way stopcock	5802.07



 OBSTETRICS NEONATOLOGY ENTERAL

 OBSTETRICS NEONATOLOGY ENTERAL

**For further information, please contact: [questions@vygon.com](mailto:questions@vygon.com)**

The specifications shown in this leaflet are for information only  
and are not, under any circumstances, of a contractual nature.

---

Vygon – 8, rue de Paris • 95440 ECOUEN • FRANCE

Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81

Export customer service: +33 (0)1.39.92.64.15

Fax.: +33 (0)1.39.92.64.44 • [www.vygon.com](http://www.vygon.com)



[www.vygon.com](http://www.vygon.com)